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### *Streptococcus pneumoniae*

#### Expert Rules and Intrinsic Resistance Tables

An MIC breakpoint of  $S \leq 0.001$  mg/L is an arbitrary, "off scale" breakpoint (corrected for mechanisms to the agent) as "Susceptible, increased exposure" (I). For these c

**MIC determination (broth microdilution according to ISO standard 20776-1)**  
**Medium:** Mueller-Hinton broth + 5% lysed horse blood and 20 mg/L  $\beta$ -NAD (MH-F br  
**Inoculum:**  $5 \times 10^5$  CFU/mL  
**Incubation:** Sealed panels, air,  $35 \pm 1^\circ\text{C}$ ,  $18 \pm 2$ h  
**Reading:** Unless otherwise stated, read MICs at the lowest concentration of the agent inhibits visible growth. See "EUCAST Reading Guide for broth microdilution" for furth  
**Quality control:** *Streptococcus pneumoniae* ATCC 49619. For agents not covered t EUCAST QC Tables.

Penicillins <sup>1</sup>	MIC breakpoints (mg/L)		
	S $\leq$	R >	ATU
Benzylpenicillin (indications other than meningitis) <sup>2</sup>	0.06	2	
Benzylpenicillin (meningitis)	0.06	0.06	

Oxacilline screening:

**Oxacillin 1  $\mu\text{g}$  zone diameter <20 mm (or benzylpenicillin MIC >0.06 mg/L)**

**Mechanism:** beta-lactam resistance detected

**Report:** resistant (R) to benzylpenicillin (meningitis) and phenoxymethylpenicillin (all indications).

For benzylpenicillin (indications other than meningitis), perform and interpret MIC according to breakpoints.

For other beta-lactam agents, see below.

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### **Pneumonia caused by *S. pneumoniae*: breakpoints are related to dosage:**


For a dose of 1.2 g (2 MU) x 4 iv, isolates with MIC  $\leq 0.5$  mg/L are susceptible.

For a dose of 2.4 (4 MU) g x 4 iv or 1.2 g (2 MU) x 6 iv, isolates with MIC  $\leq 1$  mg/L are susceptible.

For a dose of 2.4 g (4 MU) x 6 iv, isolates with MIC  $\leq 2$  mg/L are susceptible.

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**Correlation of PCG MICs for gradient tests vs. BMD  
EDL study**

Etest: BioMérieux en MTS: Liofilchem

**PCG MICs for gradient tests vs. BMD MICs  
EDL study**

**Etest, BBL MH-F**

>2 dilutions lower	0
2 dilutions lower	2
1 dilution lower	18
Identical	33
1 dilution higher	4
2 dilutions higher	0
>2 dilutions higher	0

**Etest, Oxoid MH-F**

>2 dilutions lower	1
2 dilutions lower	7
1 dilution lower	28
Identical	20
1 dilution higher	1
2 dilutions higher	
>2 dilutions higher	

**MTS, BBL MH-F**

>2 dilutions lower	
2 dilutions lower	6
1 dilution lower	33
Identical	17
1 dilution higher	1
2 dilutions higher	
>2 dilutions higher	

**MTS, Oxoid MH-F**

>2 dilutions lower	
2 dilutions lower	19
1 dilution lower	24
Identical	12
1 dilution higher	2
2 dilutions higher	
>2 dilutions higher	

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