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Algemene interpretatie

- Negative screening test:
 - MIC below or equal to or zone diameter above or equal to the susceptible breakpoint for the screening agent. No resistance mechanisms to the antimicrobial class detected.
- Positive screening test:
 - MIC above or zone diameter below the resistant breakpoint for the screening agent. Resistance mechanisms to the antimicrobial class detected.

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Staphylococcus spp.
Expert Rules and Intrinsic Resistance Tables

Macrolides, lincosamides and streptogramins	MIC breakpoints (mg/L)			Disk content (µg)	Zone diameter breakpoints (mm)		
	S ≤	R >	ATU		S ≥	R <	ATU
Azithromycin	2 ¹	2 ¹			Note ^A	Note ^A	
Clarithromycin	1 ¹	2 ¹			Note ^A	Note ^A	
Erythromycin	1	2		15	21	18	
Erythromycin (screen only)	1 ¹	1 ¹		15	21 ^A	21 ^A	

Getallen anders!!

- 1/A. Erythromycin can be used to screen for macrolide resistance in staphylococci. Isolates categorised as **screen negative** can be reported **susceptible to azithromycin, clarithromycin and roxithromycin**. Isolates categorised as **screen positive** should be **tested for susceptibility to individual agents or reported resistant**.

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De methode in het kort

Disk diffusion methodology

- **Medium:** Fastidious Anaerobe Agar (FAA) with 5% mechanically defibrinated horse blood
 - Agar depth 4.0 ± 0.5 mm
 - Plates must be dried prior to inoculation
- **Inoculum:** McF 1.0 (overnight culture from non-selective media)
- **Incubation:** Anaerobic environment
35-37°C for 16-20 h

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Methode (2)

Streaking of plates and application of disks

- For *Bacteroides* spp., remove excess fluid by turning the swab against the inside of the tube to avoid over-inoculation.
- Spread the inoculum evenly over the entire agar surface, ensuring that there are no gaps between streaks.
 - This is particularly important for *Cutibacterium acnes*, which grows with small colonies and poor contrast to the FAA media.
- Limit the number of disks on each plate to allow good growth and to avoid overlapping of zones.

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
 **EUCAST** EUROPEAN COMMITTEE ON ANTIMICROBIAL SUSCEPTIBILITY TESTING
European Society of Clinical Microbiology and Infectious Diseases

Reading guide

EUCAST disk diffusion for selected rapidly growing anaerobic bacteria on Fastidious Anaerobe Agar (FAA)

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

Anaerobic bacteria EUCAST Clinical Breakpoint Tables v. 12.0, valid from 2022-01-01

For species not listed below, see EUCAST Guidance Document on how to test and interpret results when there are no breakpoints
Expert Rules and Intrinsic Resistance Tables

MIC determination (agar dilution)
Medium: Fastidious Anaerobe Agar (FAA)
Inoculum: 10⁷ CFU/spot
Incubation: Anaerobic environment, 35-37°C, 48h
Reading: Unless otherwise stated, read MICs at the lowest concentration of the agent where a noticeable difference is seen in visible growth between the test and control plate.
Quality control: *Bacteroides fragilis* ATCC 25285 and *Clostridium perfringens* ATCC 13124.
Clostridium perfringens: For control of the inhibitor component of beta-lactam inhibitor combinations, see EUCAST QC Tables.
Clostridium perfringens DSM 25589 with a metronidazole 5 µg disk to monitor the anaerobic atmosphere.

Disk diffusion (EUCAST standardised disk diffusion method)
Medium: Fastidious Anaerobe Agar (FAA). The plates should be dried prior to inoculation (at 20-25°C overnight or at 35°C, with the lid removed, for 15 min).
Inoculum: McFarland 1.0
Incubation: Anaerobic environment, 35-37°C, 18±2h
Reading: Unless otherwise stated, read zone edges as the point showing no growth viewed from the front of the plate with the lid removed and with reflected light. See pictures below and the EUCAST Reading Guide for disk diffusion of anaerobic bacteria for further information.
Quality control: *Bacteroides fragilis* ATCC 25285 and *Clostridium perfringens* ATCC 13124. For control of the inhibitor component of beta-lactam inhibitor combination disks, see EUCAST QC Tables.
Clostridium perfringens DSM 25589 with a metronidazole 5 µg disk to monitor the anaerobic atmosphere.

Bacteroides spp.

Antimicrobial agent	MIC breakpoints (mg/L)			Disk content (µg)	Zone diameter breakpoints (mm)			Notes
	S ≤	R >	ATU		S ≤	R <	ATU	
Piperacillin-tazobactam	8 ¹	8 ¹		38-6	20	20		Notes Numbered notes relate to general comments and/or MIC breakpoints. Lettered notes relate to the disk diffusion method. 1. For susceptibility testing purposes, the concentration of tazobactam is fixed at 4 mg/L. 2A. The meropenem zone diameter breakpoint will detect all cfaA gene mediated carbapenem resistance in <i>Bacteroides fragilis</i> . Some isolates with an MIC of 1 mg/L may harbour the cfaA gene. 3/B. For information on how to use breakpoints in brackets, see https://www.eucast.org/eucastguidancedocuments/ . C. Examine zones carefully for colonies within zones. Colonies should be taken into account when reading.
Piperacillin-tazobactam, β-lactamase resistant	IE	IE			IE	IE		
Meropenem	1 ¹	1 ¹		10	28 ¹	28 ¹		
Clindamycin	(8) ¹	(4) ¹		2	(10) ^{1,2}	(10) ^{1,2}		
Metronidazole	4	4		5	25	25		

Prevotella spp.

Antimicrobial agent	MIC breakpoints (mg/L)			Disk content (µg)	Zone diameter breakpoints (mm)			Notes
	S ≤	R >	ATU		S ≤	R <	ATU	
Benzylpenicillin	0.5	0.5		1 unit	20	20		Notes Numbered notes relate to general comments and/or MIC breakpoints. Lettered notes relate to the disk diffusion method. 1. For susceptibility testing purposes, the concentration of tazobactam is fixed at 4 mg/L. A. Examine zones carefully for colonies within zones. Colonies should be taken into account when reading.
Piperacillin-tazobactam	0.5 ¹	0.5 ¹		38-6	26	26		
Meropenem	0.25	0.25		10	34	34		
Clindamycin	0.25	0.25		2	31 ¹	31 ¹		
Metronidazole	4	4		5	22	22		

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Penicillins ¹	MIC breakpoints (mg/L)			Disk content (µg)	Zone diameter breakpoints (mm)		
	S ≤	R >	ATU		S ≥	R <	ATU
Benzympenicillin (indications other than meningitis) ²	0.06	2			Note ^A	Note ^A	
Benzympenicillin (meningitis)	0.06	0.06			Note ^A	Note ^A	
Ampicillin (indications other than meningitis)	0.5	1		2	22	19	
Rifampicin	0.125	0.125		5	22	22	

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Neisseria meningitidis

Penicillins ¹	MIC breakpoints (mg/L)		
	S ≤	R >	ATU
Benzyloxyethylpenicillin (all indications)	0.25	0.25	

All indications

Cefotaxime (all indications) ¹	0.125	0.125	
---	-------	-------	--

All indications

Ceftriaxone (all indications including prophylaxis) ¹	0.125	0.125	
--	-------	-------	--

All indications including prophylaxis

Meropenem (indications other than meningitis)	Note ²	Note ²
Meropenem (all indications) ^{1,2}	0.25	0.25

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Minocycline (prophylaxis only)	1 ¹	1 ¹	
Tetracycline (screen only)	2 ¹	2 ¹	

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Oude 2021 schema Pneumokokken

Streptococcus pneumoniae EUCAST Clinical Breakpoint Tables v. 11.0, valid from 2021-01-01

Expert Rules and Intrinsic Resistance Tables

An MIC breakpoint of ≤ 0.061 mg/L is an arbitrary, "off scale" breakpoint (corresponding to a zone diameter breakpoint of " ≥ 20 mm") which categorises wild-type organisms (organisms without phenotypically detectable resistance mechanisms to the agent) as "Susceptible, increased exposure" (I). For these organism-agent combinations, never report "Susceptible, standard dosing regimen" (S).

Screening for beta-lactam resistance in *S. pneumoniae*

Oxacillin 1 µg disk test or benzylpenicillin MIC

Oxacillin zone diameter ≥ 20 mm or benzylpenicillin MIC ≤ 0.06 mg/L
Excludes all beta-lactam resistance mechanisms

Report susceptible (S) to any beta-lactam agents for which clinical breakpoints are available, including those with "Note", except for cefaclor, which if reported, should be reported "susceptible, increased exposure" (I)

Benzylpenicillin (meningitis) and phenoxymethylpenicillin (all indications)

Report resistant (R)

Benzylpenicillin (indications other than meningitis)

Determine the MIC and interpret according to the clinical breakpoints

Oxacillin zone diameter < 20 mm or benzylpenicillin MIC > 0.06 mg/L
Beta-lactam resistance mechanism detected

See the EUCAST warning on the use of benzylpenicillin gradient tests at <http://www.eucast.org/warnings/>.

Ampicillin, amoxicillin and piperacillin (without and with beta-lactamase inhibitor), cefepime, cefotaxime, ceftaroline, ceftobiprole and ceftazidime

Oxacillin zone ≥ 8 mm

Report susceptible (S)

For intravenous ampicillin, amoxicillin and piperacillin (without and with inhibitor), infer susceptibility from ampicillin

Oxacillin zone < 8 mm

Report susceptible (S)

For oral amoxicillin (without and with inhibitor), see breakpoint recommendations

Other beta-lactam agents

Determine the MIC and interpret according to the clinical breakpoints

For cefepime, cefotaxime, ceftaroline, ceftobiprole and ceftazidime, determine the MIC and interpret according to the clinical breakpoints

* In meningitis confirm by determining the MIC for the agent considered for clinical use.

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Nieuwe schema pneumokokken

Streptococcus pneumoniae EUCAST Clinical Breakpoint Tables v. 12.0, valid from 2022-01-01

Expert Rules and Intrinsic Resistance Tables

***Streptococcus pneumoniae*: Flow chart based on the oxacillin screen test for beta-lactam resistance mechanisms to reduce the number of specific tests for beta-lactam agents**

Oxacillin 1 µg zone diameter ≥ 20 mm (or benzylpenicillin MIC ≤ 0.06 mg/L)

Mechanism: excludes all beta-lactam resistance mechanisms

Report: susceptible (S) to beta-lactam agents for which clinical breakpoints are available, including those with "Note", and those with meningitis breakpoints. **Exception:** Cefaclor is reported "susceptible, increased exposure" (I).

No further testing required.

Oxacillin 1 µg zone diameter < 20 mm (or benzylpenicillin MIC > 0.06 mg/L)

Mechanism: beta-lactam resistance detected

Report: resistant (R) to benzylpenicillin (meningitis) and phenoxymethylpenicillin (all indications). For benzylpenicillin (indications other than meningitis), perform and interpret MIC according to breakpoints.

For other beta-lactam agents, see below.

Oxacillin 1 µg zone diameter 9-19 mm

Report susceptible (S) without further testing to: ampicillin, amoxicillin and piperacillin (without and with beta-lactamase inhibitor), cefepime, cefotaxime, ceftaroline, ceftobiprole, ceftazidime, imipenem and meropenem.

For other beta-lactam agents, perform susceptibility testing for the relevant agent and interpret according to breakpoints.

This guidance is also valid for meningitis breakpoints.


Oxacillin 1 µg zone diameter < 9 mm

Perform susceptibility testing for the relevant agent and interpret according to breakpoints.

This guidance is also valid for meningitis breakpoints.

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H+
MC

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Version 3.0, valid from 2022-01-18

Species	Drug	Clinical Breakpoints (mg/L)					Recommendation for area of technical uncertainty (ATU) results
		ECOFF (mg/L)	S ≤	I	R >	ATU	
<i>C. albicans</i>	Amphotericin B	1	1		1		
	Anidulafungin	0.03	0.03		0.03		
	Micafungin	0.016	0.016		0.016	0.03	If S to anidulafungin, report as S and add the following comment: Isolates susceptible to anidulafungin with micafungin MIC of 0.03 mg/L do not harbour an <i>fkx</i> mutation conferring resistance to the echinocandins. If not S to anidulafungin, report as R and refer to reference laboratory for <i>fkx</i> sequencing and confirmation of MICs.
	Fluconazole	0.5	2	4	4		
	Itraconazole	0.03	0.06		0.06		
	Posaconazole	0.06	0.06		0.06		
	Voriconazole	0.03	0.06	0.125-0.25	0.25		
<i>C. dubliniensis</i>	Amphotericin B	0.25	1		1		
	Anidulafungin	ND					
	Micafungin	ND					
	Fluconazole	(0.5) ^a	2	4	4		
	Itraconazole	0.06	0.06		0.06		
	Posaconazole	0.06	0.06		0.06		
	Voriconazole	0.03	0.06	0.125-0.25	0.25		
<i>C. glabrata</i>	Amphotericin B	1	1		1		
	Anidulafungin	0.06	0.06		0.06		
	Micafungin	0.03	0.03		0.03		
	Fluconazole	16	0.001	≤16	16		
	Itraconazole	2	ND		ND		
	Posaconazole	1	ND		ND		
	Voriconazole	1	ND		ND		
<i>C. krusei</i>	Amphotericin B	1	1		1		
	Anidulafungin	0.06	0.06		0.06		
	Micafungin	0.25	ND		ND		
	Fluconazole	128	ND		ND		
	Itraconazole	1	ND		ND		
	Posaconazole	0.5	ND		ND		
	Voriconazole	1	ND		ND		

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